

	GS1000	GS1200	GS1500	GS1800
Dimensioni esterne	1000x750x1500 mm	1200x750x1500 mm	1500x750x1500 mm	1800x750x1500 mm
Dimensioni interne	950x600x750 mm	1150x600x750 mm	1450x600x750 mm	1750x600x750 mm
Apertura lavoro max.	780x380 mm	980x380 mm	1280x380 mm	1580x380 mm
Apertura servizio max.	780x650 mm	980x650 mm	1280x650 mm	1580x650 mm
Illuminazione(220V/50Hz) fluorescente	2x18W	2x18W	2x36W	2x36w
Alimentaizone (220V/50Hz)	3A	3A	2.6A	2.6A
Portata nominale	700 m ³ /h	800 m ³ /h	950 m ³ /h	1300 m ³ /h
Numero di colonne filtranti	1	2	2	3
Peso (filtri esclusi)	140 Kg	174 Kg	200 Kg	240 Kg
Imballo	1200x800x1600 mm	1400x800x1600 mm	1600x800x1600 mm	2000x800x1600 mm