

Malt Extract Agar with Chloramphenicol

Shorthand Symbol:	MEAC
Item Number:	50-1556
Format:	250ml flask containing 200 ml medium
Colour:	Yellowish
Storage:	Constant temperature, protected from light, at 4-25 °C
Shelf-life:	?
pH:	5,6 ± 0.2

Intended use and fields of application (Specification)

Detection, isolation and determination of the germ count of yeasts and fungi in foodstuffs and other test materials.

The medium is rendered selective for yeast and fungi by the addition of Chloramphenicol.

Composition (g/l)

Malt extract	30
Chloramphenicol	0,05
Peptone	5
Agar	15

Xebios is certified according to DIN EN ISO 9001 as well as DIN EN ISO 13485. The microbiological quality control is carried out by our laboratory accredited according to DIN EN ISO / IEC 17025 (registration number of the certificate: D-PL-17997-01-00).