



Microfuge 20, 220-240 V 50-60 Hz, with FA241.5P, IVD Cod. B30140



Maximum Speed	15,000 rpm (in 100 rpm increments)
Maximum RCF	20,627 x g
Maximum Capacity	36 x 1.5/2.0 mL
Run Time	99 min, 59 sec, continuous run, short run
User-Defined Programs	10
Noise Level at Max. Speed	< 59 dBa
Accel/Decel Profiles	fast soft
Ambient Temperature Range	5° to 40°C
Maximum Heat Dissipation into Room	1433 Btu/h (420 W)
Consumo Energetico	170 (W)
Dimensions (in)	12.2W x 16.5D x 10.7H
Dimensions (cm)	31.0W x 42.0D x 27H
Weight	13 kg (28.6 lb)
Electrical Requirements	220–240 VAC, 50/60 Hz



- Intuitive user interface displays both speed (rpm) and RCG (xg)
- Decrease setup time with user-defined programs for common protocols