

	S1000	S1200	S1500	S1800
Dimensioni esterne	1000x750x1500 mm	1200x750x1500 mm	1500x750x1500 mm	1800x750x1500 mm
Dimensioni interne	950x600x1150 mm	1150x600x1150 mm	1450x600x1150 mm	1750x600x1150 mm
Apertura lavoro max.	780x380 mm	980x380 mm	1280x380 mm	1580x380 mm
Apertura servizio max.	780x750 mm	980x750 mm	1280x750 mm	1580x750 mm
Illuminazione(220V/50Hz) fluorescente	2x18W	2x18W	2x36W	2x36w
Alimentaizone (220V/50Hz)	1.5A	1.5A	1.5A	1.5A
Portata nominale	700 m ³ /h	800 m ³ /h	950 m ³ /h	1300 m ³ /h
Peso	116 Kg	132 Kg	157 Kg	182 Kg
Imballo	1200x800x1600 mm	1400x800x1600 mm	1600x800x1600 mm	2000x800x1600 mm