

Specifiche

Volume	Channel	Test volume	Increments	Inaccuracy	Imprecision	Cat. No.
0.5-10 µL	8	0.50 µL	0.1 µL	4.00 %	4.00 %	8-201-20-9
		5 µL		2.50 %	2.50 %	
		10 µL		1.50 %	1.50 %	
5-50 µL	8	5 µL	0.5 µL	3.00 %	2.00 %	8-202-20-9
		25 µL		1.50 %	1.00 %	
		50 µL		1.00 %	0.70 %	
10-100 µL	8	10 µL	0.5 µL	1,50 %	0,75 %	8-204-20-9
		50 µL		1,00 %	0,50 %	
		100 µL		1,00 %	0,25 %	
30-300 µL	8	30 µL	1.0 µL	1.50 %	0.75 %	8-203-20-9
		150 µL		1.00 %	0.50 %	
		300 µL		0.80 %	0.25 %	
0.5-10 µL	12	1 µL	0.1 µL	4.00 %	4.00 %	8-401-20-9
		5 µL		2.50 %	2.50 %	
		10 µL		1.50 %	1.50 %	
5-50 µL	12	5 µL	0.5 µL	3.00 %	2.00 %	8-402-20-9
		25 µL		1.50 %	1.00 %	
		50 µL		1.00 %	0.70 %	
10-100 µL	12	10 µL	1,00 µL	1,50 %	0,75 %	8-404-20-9
		50 µL		1,00 %	0,50 %	
		100 µL		1,00 %	0,25 %	
30-300 µL	12	30 µL	1.0 µL	1.50 %	0.75 %	8-403-20-9
		150 µL		1.00 %	0.50 %	
		300 µL		0.80 %	0.25 %	