

Technical Specification



Product: VeriGuard IPA
Non-woven
Wipe
Code: WSVG04023



Product	Pre-Saturated Wipe
Brand	VeriGuard
Material	Polycellulose
Size	230mm x 230mm
Colour	White

PRODUCT INFORMATION

Former Code(s)	vG-3-100-S 71103S
Product Description	Sterile white hydro-entangled polycellulose pouch wipe impregnated with 70% USP IPA v/v in 30% DI water Wipe size 230mm x 230mm
Intended Use	Sterile hard surface IPA wipe
Pack Format	100 wipes per pouch double-bagged in polythene packaging 8 pouches per carton
Carton Label Information and Traceability	<ul style="list-style-type: none"> Product description Product code Batch number Expiry date Pack quantity Irradiation number & indicator
Carton Specification	Dimensions (mm) = 270 x 200 x 490 Weight (kg) = 8



TECHNICAL INFORMATION

Materials and Performance	Wipe size and dry weight = 230mm x 230mm, 68±3g/m ² Wipe thickness = 0.30±0.02mm ISO 9073-6 absorption capacity and rate = 366-448% in 2.48 (±0.25) seconds Hydroentangled blend of 45% polyester and 55% cellulose 'C-folded' wipes contained in resealable pouch impregnated with 70% isopropyl alcohol and 30% DI water
Filtration	Contents 0.2µm filtered
Particle Release (≥5µm)	IEST-RP-CC004.3: Particles >0.5µm /m ² (x10 ⁶) ≤878 IEST-RP-CC004.3: Fibres >100µm/m ² ≤1.05
Packaging	Primary packaged in polythene bags in ISO 14644 Class 5 (at rest) cleanroom
Sterilisation Method	Sterilised by gamma irradiation (25kGy dose minimum)
Quality Assurance	Product supplied under BS EN ISO 9001:2008 registered systems
Storage Conditions	Store in a cool, dry, well-ventilated area, out of direct sunlight and away from sources of heat or ignition
Hazard Statement	See MSDS for details
Shelf Life	2 years from date of manufacture

Contact Details

Micronclean Skegness

Roman Bank
Skegness
Lincolnshire
PE25 1SQ
UK

Tel: +44 (0) 1754 767377
orders@micronclean.co.uk

Micronclean Newbury

Faraday Road
Newbury
West Berkshire
RG14 2AD
UK

Tel: +44 (0) 1635 234050
consumables@micronclean.biz

Micronclean Products Ltd.

Rookery House
The Guineas
Newmarket, Suffolk
CB8 8EQ
UK

Tel: +44 (0) 1842 822150
sales@micronclean.com