

# Technical Specification



**Product:** VeriGuard IPA  
Non-woven  
Wipe  
**Code:** WSVG05020



<b>Product</b>	Pre-Saturated Wipe
<b>Brand</b>	VeriGuard
<b>Material</b>	Thermobonded nonwoven PP
<b>Size</b>	200mm x 195mm
<b>Colour</b>	White

## PRODUCT INFORMATION

Former Code(s)	vG-1-200-S
Product Description	Continuous filament thermobonded nonwoven 100% PP hydrophilic IPA tub wipes Wipe size 200mm x 195mm
Intended Use	Hard surface IPA wipe
Pack Format	200 wipes per tub double-bagged 12 tubs per carton
Carton Label Information and Traceability	<ul style="list-style-type: none"> <li>Product description</li> <li>Product code</li> <li>Batch number</li> <li>Expiry date</li> <li>Pack quantity</li> <li>Irradiation number &amp; indicator</li> </ul>
Carton Specification	Dimensions (mm) = 460 x 340 x 370 Weight (kg) = 10



**TECHNICAL INFORMATION**

Materials and Performance	Thermobonded nonwoven polypropylene wipe Wipe size and dry weight = <i>195mm x 200mm, 18g/m<sup>2</sup> (±7.5%)</i> Wipe thickness = <i>0.2mm (&gt;0.15:&lt;0.3mm)</i> Absorption capacity = <i>&gt;450% average (&gt;330% minimum)</i> Impregnated with 70% isopropyl alcohol and 30% DI water Log roll contained in HDPE tub with resealable flip-top lid
Filtration	Contents 0.2µm filtered
Particle Release (≥5µm)	TEST-RP-CC003.3 Helmke Drum Test at 10rpm (0.5µm particles) = <i>41 counts/wipe/ft<sup>3</sup>/min</i>
Packaging	Individually double packaged in polythene bags
Sterilisation Method	Sterilised by gamma irradiation (25kGy dose minimum)
Quality Assurance	Product supplied under BS EN ISO 9001:2008 registered systems
Storage Conditions	Store in a cool, dry, well-ventilated area, out of direct sunlight and away from sources of heat or ignition
Hazard Statement	See MSDS for details
Shelf Life	2 years from date of manufacture

**Contact Details****Micronclean Skegness**

Roman Bank  
Skegness  
Lincolnshire  
PE25 1SQ  
UK

Tel: +44 (0) 1754 767377  
orders@micronclean.co.uk

**Micronclean Newbury**

Faraday Road  
Newbury  
West Berkshire  
RG14 2AD  
UK

Tel: +44 (0) 1635 234050  
consumables@micronclean.biz

**Micronclean Products Ltd.**

Rookery House  
The Guineas  
Newmarket, Suffolk  
CB8 8EQ  
UK

Tel: +44 (0) 1842 822150  
sales@micronclean.com