


# Technical Specification



**Product:** IPA 5 Litre  
**Code:** ZSAL64005



<b>Product</b>	Cleaning Agent
<b>Brand</b>	Compliance 100
<b>Solution</b>	IPA / WFI
<b>Size</b>	5 litres

## **PRODUCT INFORMATION**

Product Description	Sterile double-bagged 5 litre pour bottles containing IPA 70/30
Intended Use	Hard surface alcohol disinfectant suitable for ISO 14644 Class 4 cleanrooms
Pack Format	2 x 5 litre bottles per carton
Carton Label Information and Traceability	<ul style="list-style-type: none"> <li>▪ Product description</li> <li>▪ Product code &amp; barcode</li> <li>▪ Batch number</li> <li>▪ Expiry date</li> <li>▪ Size</li> <li>▪ Pack quantity</li> <li>▪ Irradiation number</li> <li>▪ Irradiation indicator</li> </ul>
Carton Specification	Dimensions (mm) = 305 x 229 x 305 Weight (kg) = 10

**TECHNICAL INFORMATION**

Materials and Performance	Natural HDPE 5 litre bottle with tamper-evident PP vented cap filled with 70% isopropyl alcohol / 30% water for injection
Filtration	Contents 0.2µm filtered
Packaging	Primary and secondary packaged in 38µm polythene bags in ISO 14644 Class 4, EU GMP Grade B cleanroom
Sterilisation Method	Validated sterilisation by gamma irradiation (min. 25 kGy)
Quality Assurance	Product manufactured under BS EN ISO 9001:2008 registered systems
Storage Conditions	Store in a cool, dry, well-ventilated area away from direct sunlight, heat and naked flames
Hazard Statement	See MSDS for detailed information
Shelf Life	2 years from date of manufacture

**Contact Details****Micronclean Skegness**

Roman Bank  
Skegness  
Lincolnshire  
PE25 1SQ  
UK

Tel: +44 (0) 1754 767377  
orders@micronclean.co.uk

**Micronclean Newbury**

Faraday Road  
Newbury  
West Berkshire  
RG14 2AD  
UK

Tel: +44 (0) 1635 234050  
consumables@micronclean.biz

**Micronclean Products Ltd.**

Rookery House  
The Guineas  
Newmarket, Suffolk  
CB8 8EQ  
UK

Tel: +44 (0) 1842 822150  
sales@micronclean.com