



## **IMS 70/30 500ml TRIGGER SPRAY SPECIFICATION**

<b>Product Code :</b>	ZSAL65499
<b>Product Description :</b>	Sterile cleanroom trigger spray containing 70% denatured ethanol / 30% water for injection
<b>Pack Size :</b>	24 x 500ml bottles per carton

### **Technical Specification**

<b>Packaging :</b>	Natural HDPE Bottle with tamper-proof trigger spray
<b>Bottle Contents :</b>	500ml
<b>Alcohol :</b>	70% USP Denatured Ethanol / 30% USP WFI
<b>Decontamination :</b>	Product filled in an ISO 14644-1 Class 5 (in use) unidirectional flow cabinet, secondary packaging in an ISO 14644-1 Class 7 (in use) cleanroom
<b>Sterilisation :</b>	Validated sterilisation by gamma irradiation (minimum 25 kGy)
<b>Shelf-Life :</b>	Two years from manufacturing date
<b>Quality Assurance :</b>	Product manufactured under BS EN ISO 9001:2008 registered systems

**Issue no 1 30/06/10**