



COMPLIANCE 100® ALPHA 5-LITRE RFU

TECHNICAL / PRODUCT SPECIFICATION

Micronclean Product Code	ZSAP03005
Product Description	Sterile 5-litre Compliance 100® Alpha pour bottles
Bottle Description	Natural HDPE bottle with tamper-evident PP vented cap
Packaging Configuration	Individually double-bagged in polythene packaging
Pack / Carton Label Information	<ul style="list-style-type: none"> • Product description • Micronclean product code and barcode • Irradiation number • Expiry date • Pack quantity • Irradiation indicator
Disinfectant Specification	Compliance 100® Alpha broad spectrum disinfectant
Decontamination	Product filled and secondary packaged in an EU GMP Grade B and ISO 14644-1 Class 6 (at rest) cleanroom
Filtration	Filtered to 0.2µm
Carton Configuration	2 x 5 litre bottles per carton
Sterilisation Method	Gamma irradiation
Sterilisation Sub-contractor	Synergy Healthcare
Shelf Life	2 years
Storage	Store in a cool, dry place away from direct sunlight, heat and naked flames