

Technical Specification



Product: GuardMark 2
Sterile
Marker Pen

Codes:
ZSBK03010 (black)



Product	Stationery
Brand	GuardMark 2
Material	Cleanroom Marker Pen
Colour	Black

PRODUCT INFORMATION

Product Description

Sterile individually packed marker pens

Intended Use

Sterile cleanroom marker pens for marking on plastic, glass and metal

Pack Format

10 individually double-packaged marker pens per outer packaging

Carton Label Information and
Traceability

- Product description
- Colour
- Pack quantity
- Product code & barcode
- Batch number
- Expiry date
- Pack quantity
- Irradiation indicator

Carton Specification

Dimensions (mm) = 190 x 101 x 51mm
Weight (kg) = 0.16



TECHNICAL INFORMATION

Materials and Performance

Materials consist of

- Barrel, Cap and Plug – Polypropylene
- Reservoir and nib – Polyester
- Lab Marker Ink – see attached MSDS

Size = 134mm x 17mm

Used to make permanent markings on plastic, metal and porcelain lab ware. The unique ink has the ability to write on lab ware even when the surface is damp. These lab markers have been tested against water and isopropyl alcohol and have shown resistance to both.

Packaging

Primary packaged in polythene bags and secondary packaged in cardboard carton in ISO 14644 Class 6 (at rest) cleanroom

Sterilisation Method

Sterilised by gamma irradiation (25kGy dose minimum)

Quality Assurance

Product supplied under BS EN ISO 9001:2008 registered systems

Storage Conditions

Store in ambient conditions.

Hazard Statement

Non-hazardous

Shelf Life

2 years from date of manufacture

Contact Details

Micronclean Skegness

Roman Bank
Skegness
Lincolnshire
PE25 1AQ
UK

Tel: +44 (0) 1754 767377
orders@micronclean.co.uk

Micronclean Products Ltd.

Rookery House
The Guineas
Newmarket, Suffolk
CB8 8EQ
UK

Tel: +44 (0) 1842 822150
sales@micronclean.com